


Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.844 (0.683 to 1.005)	0.725 (0.569 to 0.881)
Intercept:	0.01 (-0.01 to 0.02)	0.02 (0.00 to 0.03)
Std Err Est:	0.04	0.04

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.8361	Points (Plotted/Total): 40/40
Bias: 0.00	Outliers: Not tested
XMean ± SD: 0.07 ± 0.08	Degrees Freedom: 38
YMean ± SD: 0.06 ± 0.07	Scatter Plot Bounds: None
Std Dev Diff: 0.04	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	14 Oct 2002	14 Oct 2002
Rep SD:	1	1
Result Ranges:	0.0 to 0.4	0.0 to 0.3
Units:	x10 ⁹ /L	10 ⁹ /L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	0.0	0.0	0.0	S00015	0.0	0.0	0.0	S00029	0.1	0.1	0.0
S00002	0.1	0.1	0.0	S00016	0.1	0.1	0.0	S00030	0.0	0.0	0.0
S00003	0.0	0.0	0.0	S00017	0.1	0.1	0.0	S00031	0.1	0.1	0.0
S00004	0.1	0.1	0.0	S00018	0.1	0.1	0.0	S00032	0.2	0.2	0.0
S00005	0.1	0.0	-0.1	S00019	0.1	0.1	0.0	S00033	0.0	0.0	0.0
S00006	0.0	0.0	0.0	S00020	0.0	0.1	0.1	S00034	0.0	0.0	0.0
S00007	0.1	0.1	0.0	S00021	0.1	0.1	0.0	S00035	0.1	0.1	0.0
S00008	0.0	0.1	0.1	S00022	0.0	0.0	0.0	S00036	0.4	0.3	-0.1
S00009	0.1	0.1	0.0	S00023	0.0	0.1	0.1	S00037	0.0	0.0	0.0
S00010	0.1	0.0	-0.1	S00024	0.0	0.0	0.0	S00038	0.1	0.0	-0.1
S00011	0.1	0.1	0.0	S00025	0.0	0.0	0.0	S00039	0.0	0.0	0.0
S00012	0.1	0.1	0.0	S00026	0.0	0.0	0.0	S00040	0.1	0.1	0.0
S00013	0.1	0.1	0.0	S00027	0.1	0.1	0.0				
S00014	0.0	0.0	0.0	S00028	0.0	0.0	0.0				

Values with an "X" were excluded from the calculations.