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## Test strips for rapid determination of blood, protein, ascorbic acid, glucose and pH-value in urine

### Use

Screening test for detection of diabetes and diseases of kidney and urinary tract.

**Only for use by qualified personnel.**

### Instructions for use

Dip the test strip for approximately 1 second into the fresh urine. Draw it across the rim of the container to remove excess urine. After 30 to 60 seconds compare the test strip with the color scale. The best time for comparison is after 30 seconds. Color changes that take place after more than 2 minutes are of no significance. When tested the urine should not be older than 2 hours.

### Principle

**Blood:** The detection is based on the pseudoperoxidative activity of hemoglobin and myoglobin, which catalyze the oxidation of an indicator by an organic hydroperoxide producing a green color.

**Protein:** The test is based on the “protein error” principle of indicators. The test zone is buffered to a constant pH value and changes color from yellow to greenish blue in the presence of albumin. Other proteins are indicated with less sensitivity.

**Ascorbic acid:** The detection is based on the decoloration of Tillmans reagent. In the presence of ascorbic acid a color change takes place from blue to red.

**Glucose:** The detection is based on the glucoseoxidase-peroxidase-chromogen reaction. Apart from glucose, no other compound in urine is known to give a positive reaction.

**pH:** The test paper contains indicators which clearly change color between pH 5 and pH 9 (from orange to green to turquoise).

### Evaluation – Sources of Error

**Blood:** The minimum sensitivity of the test strip is 5 erythrocytes/µL urine corresponding to approx. 0.015 mg hemoglobin/dL urine. Intact erythrocytes are indicated by flecky discolorations of the test field. The color fields correspond to the following values:

0 (negative), ca. 5–10, ca. 50, ca. 250 Ery/µL resp.

hemoglobin concentration out of ca. 10, ca. 50, ca. 250 Ery/µL

Normal concentrations of ascorbic acid (< 40 mg/dL) do not influence the test results. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.

**Protein:** The minimum sensitivity of the test strip is 10 mg protein/dL urine. The color fields correspond to the following ranges of albumin concentrations:

negative, 30, 100, 500 mg/dL or

negative, 0,3, 1,0, 5,0 g/L

Falsely positive results are possible in alkaline urine samples (pH > 9), after infusions with polyvinylpyrrolidone (blood substitute), after intake of medicaments containing quinine and also by disinfectant residues in the urine sampling vessel. The protein coloration may be masked by the presence of medical dyes (e.g. methylene blue) or beetroot pigments.

**Ascorbic acid:** The color fields correspond to the following values:

0 (negative), 10(+), 20(++) mg/dL or

0 (negative), 0,6(+), 1,1(++) mmol/L

Only for information!

**Glucose:** Pathological glucose concentrations are indicated by a color change from green to bluish green. Yellow or greenish test fields should be considered negative or normal. The color fields correspond to the following ranges of glucose concentrations:

neg. (yellow), neg. or normal (greenish), 50, 150, 500, ≥ 1000 mg/dL or

neg. (yellow), neg. or normal (greenish), 2,8, 8,3, 27,8, ≥ 55,5 mmol/L

The influence of ascorbic acid (vitamin C) has been largely eliminated. An inhibitory effect is produced by gentisic acid. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.

**pH:** The pH value of fresh urine of healthy people varies between pH 5 and pH 6. The color scale gives a clear distinction of pH value between pH 5 and pH 9.

### Quality Control in professional use

The performance of the test strips should be confirmed by use of positive and negative control solutions. Positive and negative controls should be analyzed once a day, whenever a new bottle of strips is opened, whenever a new lot of strips is started, and every 30 days to check storage conditions. Each laboratory should establish its own goals for adequate standards of performance, and should question handling and testing procedures if these standards are not met.

### Reactive ingredients

(Quantity resp. activity/cm<sup>2</sup> at time of impregnation)

<b>Blood:</b>		<b>Glucose:</b>	
tetramethylbenzidine	31 µg	glucose oxidase	7 U
cumene hydroperoxide	315 µg	peroxidase	1 U
		tetramethylbenzidine	96 µg
<b>Protein:</b>		<b>pH:</b>	
tetrabromophenol blue	10 µg	methyl red	3 µg
		bromothymol blue	10 µg
<b>Ascorbic acid:</b>			
2,6-dichlorophenolindophenol	7 µg		

### Directions

In any case, in order to establish a final diagnosis and prescribe an appropriate therapy, the results obtained with test strips should be verified with other medical results.

The effect of medicaments or their metabolic products on the test is not known in all cases. In case of doubt it is recommended not to take the medicaments and then repeat the test.

Only use well washed and clean vessels for urine collection. The presence of usual urine preservatives will not affect the test results.

Remove only as many test strips as are required, and reseal the container immediately after use. Do not touch the test pads. Avoid exposing the strips to sunlight and moisture. Store the container below + 30 °C in a dry place. The test strips are stable, when stored properly up to the date of expiry indicated.

The caps contain a non-poisonous and harmless desiccant. In case this desiccant is swallowed accidentally, then drink plenty of water.

Explanation of symbols can be found in the package insert.

Disposal: Please dispose all used test strips in accordance with your local laws and regulations.

Package units: Tubes of 50 and 100 test strips

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