

TECHNICAL DATA SHEET

DeltaSwab tube with AMIES media and flocked swab

code 304287

code	description	break point	case quantity	case weight (kg)	case volume (m ³)	cases pallet
304287	2 ml Amies + Standard flocked swab	80 mm	6x100	5.70	0.056	32



Tube with media: Directive 98/79/CE. "In vitro" diagnostic medical devices
 Swab: Directive 93/42/CEE. Medical devices

DESCRIPTION AND SPECIFIC USE

The kit includes:

- Tube with Amies medium
- Flocked swab
- Peel-pack package

Swab: sterile by radiation (R).

Tube with media: sterile by radiation (R).

Particularly suited for microbiological sample collection and transport. It maintains the reliability of Aerobes, Facultative Anaerobes and Anaerobes Bacteria for a minimum of 48 hours and 24 hours with Fastidious Bacteria.

It is designed for traditional (manual) sample streaking and it is also compatible with automated sample streaking devices of the market. It is also appropriate in molecular diagnostic techniques and for direct extension on glass slides.

MEASUREMENTS

	length	diameter
Tube	82.6 mm	13.0 mm
Cap	-	16.6 mm
Swab	152.0 mm	-

WRITTEN


David Campos
Marketing Department

VERIFIED


Anna Mir
Technical, Development and Quality Director

AUTHORISED


Mònica Torras
Sales & Marketing Director

TECHNICAL DATA SHEET

DeltaSwab tube with AMIES media and flocked swab

code 304287


FEATURES OF DELTASWAB SYSTEM

- Compatible with the new inoculation automation and sample streaking equipment.
- Compatible with molecular diagnostic techniques.
- Facilitates the collection, transport and subsequent treatment of microbiological samples.
- Greater reliability for sample recovery: high capacity for the absorption and elution of the sample from the new flocked swab Deltaswab.
- Greater feasibility of the sample as it is totally suspended in the medium.
- Allows different homogeneous inocula to be obtained from the same sample.
- It adapts to any working protocol.
- Minimises any possible cross-contamination: less handling and maximum water tightness.
- More comfort for the user: conical base which facilitates stirring and greater stability of the tube with skirt base.
- More comfort for the patient, due to the softer swab covering.
- Standardises the receipt of samples at the microbiology labs.
- Facilitates manual streaking.
- Allows direct extension on the glass slide as the medium does not contain agar.
- Ensures proper sample transport and storage both at room temperatures (20 °C – 25 °C) and at refrigeration temperatures (4 °C – 8 °C).

FEATURES OF FLOCKED SWABS

- Maximum reliability and convenience of the breakpoint.
- Support made of polystyrene --> Optimum safety.
- Head swab made of Flocked polyester --> increased absorption and elution of the sample on the media.
- Maximum adaptability: wide range
- Maximum comfort: anatomical and ergonomic design.

INSTRUCTIONS



1. Open peel-pack using both hands pulling opposite sides.
2. Collect the sample with the swab.
3. Open the tube with the other hand and place the swab inside so that it is covered with the media.
4. Align the breakpoint of the swab with the top of the tube lightly pressing the swab downwards. Break the swab by the breakpoint by supporting it on the inner edge of the tube.
5. Discard leftover piece of stick, screw cap tightly and shake the sample to elute it into the media.
6. Identify the sample and send it to the lab.

WRITTEN

David Campos
Marketing Department

VERIFIED

Anna Mir
Technical, Development and Quality Director


AUTHORISED

Mònica Torras
Sales & Marketing Director

TECHNICAL DATA SHEET

DeltaSwab tube with AMIES media and flocked swab

code 304287






Brand name → **δswab**

CE Classification → **Tubo con medio de transporte**
 Tube with transport media
IVD
 Hisopo recolector de muestras
 Swab, specimen collection
 Producto Sanitario - Medical Device
MDD CE 0318

← Instructions

Manufacturer → **DELTALAB, S.L.**
 Plaza de la Verneda, 1
 08191 Rubí (SPAIN)

Single use product →    → To maintain between 2 °C and 30 °C

Method of sterilisation → **STERILE IR**

Do not use if the packaging is damaged

WRITTEN	VERIFIED	AUTHORISED
 David Campos Marketing Department	 Anna Mir Technical, Development and Quality Director	 Mònica Torras Sales & Marketing Director