


## Alternate Method Comparison

X Method: BD

Y Method: Greiner



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	0.821 (0.668 to 0.975)	0.717 (0.568 to 0.867)
<b>Intercept:</b>	1.142 (0.059 to 2.226)	1.870 (0.815 to 2.926)
<b>Std Err Est:</b>	0.464	0.452

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R):	0.8477	Points (Plotted/Total):	39/39
Bias:	-0.107	Outliers:	None
XMean ± SD:	6.993 ± 0.993	Degrees Freedom:	37
YMean ± SD:	6.885 ± 0.840	Scatter Plot Bounds:	None
Std Dev Diff:	0.527		
SubRange Bounds:	None		

### Experiment Description

	X Method	Y Method
ExptDate:	11 Feb 2003	11 Feb 2003
Rep SD:	1	1
Result Ranges:	4.79 to 9.06	5.01 to 8.66
Units:	pmol/L	pmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_  
Signature

\_\_\_\_\_ Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

## Alternate Method Comparison

**X Method: BD**

**Y Method: Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	5.11	6.22	1.11	S00014	9.06	8.36	-0.70	S00027	8.14	7.99	-0.15
S00002	7.02	7.37	0.35	S00015	7.51	7.63	0.12	S00028	7.42	7.45	0.03
S00003	6.87	7.42	0.55	S00016	8.88	8.66	-0.22	S00029	8.05	7.05	-1.00
S00004	4.79	5.78	0.99	S00017	5.36	5.01	-0.35	S00030	7.42	7.42	0.00
S00005	6.51	6.87	0.36	S00018	7.53	7.31	-0.22	S00031	6.13	5.85	-0.28
S00006	6.07	6.67	0.60	S00019	5.45	5.25	-0.20	S00032	7.31	6.82	-0.49
S00007	6.59	7.79	1.20	S00020	7.11	6.64	-0.47	S00033	7.22	7.33	0.11
S00008	6.60	6.68	0.08	S00021	7.53	7.30	-0.23	S00034	7.85	7.74	-0.11
S00009	5.87	6.05	0.18	S00022	8.16	7.33	-0.83	S00035	6.62	5.78	-0.84
S00010	6.59	6.59	0.00	S00023	8.08	7.68	-0.40	S00036	7.23	6.60	-0.63
S00011	6.90	7.03	0.13	S00024	6.79	6.84	0.05	S00037	6.99	7.11	0.12
S00012	6.71	6.34	-0.37	S00025	5.31	5.21	-0.10	S00038	7.27	6.51	-0.76
S00013	6.64	6.08	-0.56	S00026	7.8	7.46	-0.34	S00039	8.22	7.3	-0.92

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.