


## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	0.964 (0.861 to 1.067)	0.915 (0.813 to 1.017)
<b>Intercept:</b>	0.464 (-0.999 to 1.928)	1.153 (-0.293 to 2.598)
<b>Std Err Est:</b>	0.804	0.794

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R): 0.9470	Points (Plotted/Total): 40/40
Bias: -0.038	Outliers: None
XMean ± SD: 13.967 ± 2.526	Degrees Freedom: 38
YMean ± SD: 13.929 ± 2.440	Scatter Plot Bounds: None
Std Dev Diff: 0.813	
SubRange Bounds: None	

### Experiment Description

	X Method	Y Method
ExptDate:	11 Feb 2003	11 Feb 2003
Rep SD:	1	1
Result Ranges:	9.65 to 21.62	9.78 to 21.75
Units:	pmol/L	pmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_

Signature

Date

## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	15.44	14.80	-0.64	S00015	14.8	14.03	-0.77	S00029	12.74	11.97	-0.77
S00002	17.12	14.80	-2.32	S00016	14.54	14.93	0.39	S00030	14.80	14.41	-0.39
S00003	21.62	21.75	0.13	S00017	10.55	11.33	0.78	S00031	19.69	19.95	0.26
S00004	15.96	13.51	-2.45	S00018	12.10	13.13	1.03	S00032	11.84	12.23	0.39
S00005	13.77	13.26	-0.51	S00019	11.71	12.36	0.65	S00033	13.26	13.13	-0.13
S00006	11.97	10.81	-1.16	S00020	12.10	12.23	0.13	S00034	12.36	11.58	-0.78
S00007	18.53	18.92	0.39	S00021	15.57	16.34	0.77	S00035	14.29	13.26	-1.03
S00008	13.00	13.51	0.51	S00022	13.00	12.87	-0.13	S00036	13.51	14.03	0.52
S00009	12.87	12.87	0.00	S00023	12.61	14.03	1.42	S00037	15.83	15.96	0.13
S00010	10.68	11.07	0.39	S00024	13.64	13.64	0.00	S00038	17.63	17.12	-0.51
S00011	13.26	13.51	0.25	S00025	9.65	9.78	0.13	S00039	12.74	12.61	-0.13
S00012	10.42	11.84	1.42	S00026	15.19	15.06	-0.13	S00040	13.90	14.41	0.51
S00013	16.6	16.86	0.26	S00027	12.23	11.58	-0.65				
S00014	12.74	13.26	0.52	S00028	14.41	14.41	0.00				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.