


Alternate Method Comparison


X Method: BD

Y Method: Greiner


Scatter Plot



Bias



Percent Bias



Regression Analysis

	Deming	Regular
Slope:	1.008 (0.976 to 1.040)	1.003 (0.971 to 1.035)
Intercept:	0.5 (-0.6 to 1.5)	0.6 (-0.5 to 1.6)
Std Err Est:	1.9	1.9

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9953	Points (Plotted/Total): 40/40
Bias: 0.7	Outliers: None
XMean ± SD: 27.1 ± 19.1	Degrees Freedom: 38
YMean ± SD: 27.7 ± 19.3	Scatter Plot Bounds: None
Std Dev Diff: 1.9	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	10 Feb 2003	10 Feb 2003
Rep SD:	1	1
Result Ranges:	5 to 111	7 to 111
Units:	U/L	U/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	17	17	0	S00015	111	111	0	S00029	23	24	1
S00002	22	22	0	S00016	31	26	-5	S00030	20	21	1
S00003	14	15	1	S00017	12	12	0	S00031	59	60	1
S00004	39	39	0	S00018	10	12	2	S00032	38	41	3
S00005	60	62	2	S00019	12	12	0	S00033	25	27	2
S00006	55	55	0	S00020	30	31	1	S00034	15	17	2
S00007	14	14	0	S00021	29	29	0	S00035	39	43	4
S00008	35	37	2	S00022	24	24	0	S00036	23	29	6
S00009	12	10	-2	S00023	25	25	0	S00037	26	28	2
S00010	35	33	-2	S00024	14	13	-1	S00038	29	30	1
S00011	12	12	0	S00025	28	26	-2	S00039	14	15	1
S00012	12	12	0	S00026	7	7	0	S00040	19	20	1
S00013	15	16	1	S00027	36	40	4				
S00014	5	7	2	S00028	36	35	-1				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.