


## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	0.977 (0.957 to 0.997)	0.975 (0.955 to 0.995)
<b>Intercept:</b>	0.12 (-0.01 to 0.25)	0.13 (0.00 to 0.27)
<b>Std Err Est:</b>	0.19	0.19

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R): 0.9980	Points (Plotted/Total): 40/40
Bias: -0.01	Outliers: None
XMean ± SD: 5.76 ± 3.06	Degrees Freedom: 38
YMean ± SD: 5.76 ± 2.99	Scatter Plot Bounds: None
Std Dev Diff: 0.20	
SubRange Bounds: None	

### Experiment Description

	X Method	Y Method
ExptDate:	10 Feb 2003	10 Feb 2003
Rep SD:	1	1
Result Ranges:	3.1 to 20.3	3.2 to 19.9
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_  
Signature

\_\_\_\_\_ Date

## Alternate Method Comparison

**X Method: BD**

**Y Method: Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	4.9	5.1	0.2	S00015	3.1	3.2	0.1	S00029	5.3	5.3	0.0
S00002	11.8	11.7	-0.1	S00016	4.7	4.8	0.1	S00030	4.6	4.4	-0.2
S00003	4.0	3.6	-0.4	S00017	5.1	5.1	0.0	S00031	4.7	4.8	0.1
S00004	4.3	4.4	0.1	S00018	4.9	4.7	-0.2	S00032	5.1	5.1	0.0
S00005	5.0	5.0	0.0	S00019	5.5	5.4	-0.1	S00033	4.5	4.5	0.0
S00006	6.6	6.5	-0.1	S00020	20.3	19.9	-0.4	S00034	4.5	4.5	0.0
S00007	5.5	5.4	-0.1	S00021	14.5	14.3	-0.2	S00035	5.0	4.7	-0.3
S00008	4.2	4.4	0.2	S00022	5.5	5.3	-0.2	S00036	6.4	6.4	0.0
S00009	4.1	3.9	-0.2	S00023	5.5	5.4	-0.1	S00037	5.5	5.4	-0.1
S00010	4.8	4.9	0.1	S00024	5.1	5.1	0.0	S00038	4.5	4.4	-0.1
S00011	4.0	4.1	0.1	S00025	5.5	5.5	0.0	S00039	4.2	4.6	0.4
S00012	4.8	4.9	0.1	S00026	4.9	5.6	0.7	S00040	5.3	5.3	0.0
S00013	6.3	6.4	0.1	S00027	5.3	5.6	0.3				
S00014	5.1	5.0	-0.1	S00028	5.7	5.7	0.0				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.