




Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.998 (0.960 to 1.036)	0.992 (0.954 to 1.029)
Intercept:	0.78 (-0.75 to 2.30)	1.04 (-0.48 to 2.56)
Std Err Est:	0.66	0.66

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R):	0.9934	Points (Plotted/Total):	40/40
Bias:	0.70	Outliers:	None
XMean ± SD:	40.07 ± 5.67	Degrees Freedom:	38
YMean ± SD:	40.77 ± 5.66	Scatter Plot Bounds:	None
Std Dev Diff:	0.65		
SubRange Bounds:	None		

Experiment Description

	X Method	Y Method
ExptDate:	14 Oct 2002	14 Oct 2002
Rep SD:	1	1
Result Ranges:	24.0 to 57.5	24.0 to 58.4
Units:	%	%
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	42.3	42.0	-0.3	S00015	38.2	39.2	1.0	S00029	41.8	41.9	0.1
S00002	35.8	37.0	1.2	S00016	46.5	47.1	0.6	S00030	40.4	40.8	0.4
S00003	44.4	44.6	0.2	S00017	38.2	38.6	0.4	S00031	39.9	40.6	0.7
S00004	42.7	43.6	0.9	S00018	42.5	43.2	0.7	S00032	57.5	58.4	0.9
S00005	40.7	41.2	0.5	S00019	38.1	39.1	1.0	S00033	39.0	40.1	1.1
S00006	43.9	45.4	1.5	S00020	37.6	38.5	0.9	S00034	39.4	39.7	0.3
S00007	42.2	43.0	0.8	S00021	43.8	45.4	1.6	S00035	40.5	40.5	0.0
S00008	47.8	48.3	0.5	S00022	45.9	46.5	0.6	S00036	24.0	24.0	0.0
S00009	39.9	41.8	1.9	S00023	46.6	45.3	-1.3	S00037	41.4	41.9	0.5
S00010	37.9	38.7	0.8	S00024	26.0	26.8	0.8	S00038	38.8	40.5	1.7
S00011	35.7	36.5	0.8	S00025	39.6	40.4	0.8	S00039	37.9	38.9	1.0
S00012	42.7	43.1	0.4	S00026	32.5	33.6	1.1	S00040	39.3	38.5	-0.8
S00013	41.5	43.2	1.7	S00027	38.8	39.1	0.3				
S00014	29.9	30.9	1.0	S00028	41.0	42.7	1.7				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.