




Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	1.007 (0.968 to 1.047)	1.000 (0.961 to 1.039)
Intercept:	0.00 (-0.53 to 0.52)	0.09 (-0.43 to 0.61)
Std Err Est:	0.23	0.23

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9930	Points (Plotted/Total): 40/40
Bias: 0.10	Outliers: None
XMean ± SD: 13.21 ± 1.88	Degrees Freedom: 38
YMean ± SD: 13.31 ± 1.89	Scatter Plot Bounds: None
Std Dev Diff: 0.22	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	14 Oct 2002	14 Oct 2002
Rep SD:	1	1
Result Ranges:	7.7 to 17.1	7.7 to 17.2
Units:	g/dL	g/dL
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Alternate Method Comparison

X Method: BD

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	13.7	13.8	0.1	S00015	12.9	13.3	0.4	S00029	14.0	14.2	0.2
S00002	11.9	12.0	0.1	S00016	15.8	16.0	0.2	S00030	12.6	12.8	0.2
S00003	14.2	14.2	0.0	S00017	11.9	12.0	0.1	S00031	12.4	12.5	0.1
S00004	14.2	14.1	-0.1	S00018	14.1	14.2	0.1	S00032	17.1	17.2	0.1
S00005	13.4	13.5	0.1	S00019	12.7	12.8	0.1	S00033	12.9	13.0	0.1
S00006	15.2	15.0	-0.2	S00020	12.8	12.9	0.1	S00034	12.9	13.0	0.1
S00007	14.3	14.4	0.1	S00021	14.9	15.4	0.5	S00035	13.5	13.5	0.0
S00008	15.9	16.2	0.3	S00022	15.3	15.6	0.3	S00036	7.7	7.7	0.0
S00009	13.2	13.3	0.1	S00023	15.9	15.0	-0.9	S00037	13.9	13.9	0.0
S00010	12.7	12.8	0.1	S00024	8.5	8.4	-0.1	S00038	13.1	13.7	0.6
S00011	11.6	11.6	0.0	S00025	13.0	13.1	0.1	S00039	12.7	12.8	0.1
S00012	14.5	14.6	0.1	S00026	10.7	10.7	0.0	S00040	12.9	12.8	-0.1
S00013	13.8	14.0	0.2	S00027	12.7	13.0	0.3				
S00014	9.4	9.5	0.1	S00028	13.7	13.9	0.2				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.