




Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	1.062 (0.905 to 1.219)	0.951 (0.798 to 1.103)
Intercept:	-7.7 (-34.5 to 19.1)	11.0 (-15.1 to 37.1)
Std Err Est:	13.5	13.1

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9006	Points (Plotted/Total): 39/40
Bias: 2.7	Outliers: None
XMean ± SD: 168.4 ± 28.2	Degrees Freedom: 37
YMean ± SD: 171.1 ± 29.8	Scatter Plot Bounds: None
Std Dev Diff: 13.0	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	11 Feb 2003	11 Feb 2003
Rep SD:	1	1
Result Ranges:	119 to 251	117 to 259
Units:	U/L	U/L
Analyst:	DIS	DIS
Comment:	SP #16 HEMOLYZED	

Accepted by: _____

Signature

Date

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias	
S00001	134	150	16	S00015		238	259	21	S00029	180	171	-9
S00002	119	117	-2	S00016	X	190	382		S00030	190	179	-11
S00003	189	199	10	S00017		167	204	37	S00031	192	186	-6
S00004	161	160	-1	S00018		134	126	-8	S00032	192	198	6
S00005	167	170	3	S00019		168	171	3	S00033	146	137	-9
S00006	150	171	21	S00020		193	207	14	S00034	152	144	-8
S00007	161	176	15	S00021		153	164	11	S00035	158	138	-20
S00008	130	149	19	S00022		144	164	20	S00036	184	173	-11
S00009	137	144	7	S00023		161	172	11	S00037	204	201	-3
S00010	191	195	4	S00024		150	171	21	S00038	166	149	-17
S00011	164	176	12	S00025		251	246	-5	S00039	174	167	-7
S00012	136	139	3	S00026		199	185	-14	S00040	162	147	-15
S00013	132	131	-1	S00027		186	184	-2				
S00014	184	194	10	S00028		167	157	-10				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.