




## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	1.276 (1.142 to 1.410)	1.195 (1.064 to 1.326)
<b>Intercept:</b>	-0.12 (-0.18 to -0.05)	-0.08 (-0.14 to -0.02)
<b>Std Err Est:</b>	0.08	0.08

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R): 0.9485	Points (Plotted/Total): 40/40
Bias: 0.00	Outliers: None
XMean ± SD: 0.43 ± 0.19	Degrees Freedom: 38
YMean ± SD: 0.44 ± 0.24	Scatter Plot Bounds: None
Std Dev Diff: 0.09	
SubRange Bounds: None	

### Experiment Description

	X Method	Y Method
ExptDate:	14 Oct 2002	14 Oct 2002
Rep SD:	1	1
Result Ranges:	0.2 to 1.2	0.2 to 1.4
Units:	10^9/L	10^9/L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	0.4	0.3	-0.1	S00015	0.3	0.3	0.0	S00029	0.7	0.8	0.1
S00002	0.8	0.7	-0.1	S00016	0.4	0.4	0.0	S00030	0.4	0.5	0.1
S00003	0.3	0.4	0.1	S00017	0.4	0.4	0.0	S00031	0.4	0.3	-0.1
S00004	0.5	0.5	0.0	S00018	0.6	0.7	0.1	S00032	1.2	1.4	0.2
S00005	0.3	0.2	-0.1	S00019	0.3	0.2	-0.1	S00033	0.3	0.3	0.0
S00006	0.3	0.3	0.0	S00020	0.4	0.3	-0.1	S00034	0.4	0.4	0.0
S00007	0.2	0.2	0.0	S00021	0.4	0.4	0.0	S00035	0.3	0.2	-0.1
S00008	0.4	0.4	0.0	S00022	0.3	0.2	-0.1	S00036	0.9	1.1	0.2
S00009	0.5	0.5	0.0	S00023	0.4	0.4	0.0	S00037	0.2	0.3	0.1
S00010	0.4	0.5	0.1	S00024	0.3	0.3	0.0	S00038	0.4	0.5	0.1
S00011	0.3	0.4	0.1	S00025	0.5	0.4	-0.1	S00039	0.3	0.3	0.0
S00012	0.6	0.6	0.0	S00026	0.4	0.3	-0.1	S00040	0.6	0.6	0.0
S00013	0.5	0.6	0.1	S00027	0.4	0.4	0.0				
S00014	0.3	0.2	-0.1	S00028	0.3	0.2	-0.1				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.