


Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.963 (0.867 to 1.058)	0.920 (0.825 to 1.015)
Intercept:	0.51 (-0.38 to 1.39)	0.90 (0.02 to 1.77)
Std Err Est:	0.31	0.31

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R):	0.9542	Points (Plotted/Total):	40/40
Bias:	0.16	Outliers:	None
XMean ± SD:	9.20 ± 1.06	Degrees Freedom:	38
YMean ± SD:	9.36 ± 1.02	Scatter Plot Bounds:	None
Std Dev Diff:	0.32		
SubRange Bounds:	None		

Experiment Description

	X Method	Y Method
ExptDate:	20 Dec 2002	20 Dec 2002
Rep SD:	1	1
Result Ranges:	5.3 to 10.8	5.3 to 10.7
Units:	fL	fL
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	10.8	10.5	-0.3	S00015	9.0	9.0	0.0	S00029	9.8	10.1	0.3
S00002	8.8	9.0	0.2	S00016	10.2	10.5	0.3	S00030	7.6	8.4	0.8
S00003	9.0	8.4	-0.6	S00017	9.2	9.7	0.5	S00031	9.3	9.9	0.6
S00004	9.9	9.7	-0.2	S00018	10.8	10.6	-0.2	S00032	8.3	9.2	0.9
S00005	9.1	9.1	0.0	S00019	8.9	9.0	0.1	S00033	9.2	9.8	0.6
S00006	8.9	9.1	0.2	S00020	9.6	9.7	0.1	S00034	8.4	8.9	0.5
S00007	10.0	9.9	-0.1	S00021	8.4	8.2	-0.2	S00035	10.1	10.1	0.0
S00008	9.6	9.7	0.1	S00022	9.2	9.0	-0.2	S00036	5.3	5.3	0.0
S00009	9.3	9.5	0.2	S00023	9.5	9.4	-0.1	S00037	8.8	8.9	0.1
S00010	8.8	9.0	0.2	S00024	7.3	7.5	0.2	S00038	8.6	8.9	0.3
S00011	10.1	10.3	0.2	S00025	8.2	8.7	0.5	S00039	10.6	10.7	0.1
S00012	10.3	10.3	0.0	S00026	9.5	9.8	0.3	S00040	10.4	10.6	0.2
S00013	10.7	10.5	-0.2	S00027	9.6	10.4	0.8				
S00014	8.1	8.2	0.1	S00028	8.8	9.0	0.2				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.