




Alternate Method Comparison

X Method: BD

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	1.075 (1.032 to 1.119)	1.067 (1.023 to 1.110)
Intercept:	-0.21 (-0.46 to 0.04)	-0.17 (-0.42 to 0.08)
Std Err Est:	0.39	0.39

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9924	Points (Plotted/Total): 40/40
Bias: 0.16	Outliers: None
XMean ± SD: 4.94 ± 2.91	Degrees Freedom: 38
YMean ± SD: 5.10 ± 3.12	Scatter Plot Bounds: None
Std Dev Diff: 0.43	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	08 Jan 2003	08 Jan 2003
Rep SD:	1	1
Result Ranges:	2.3 to 16.7	2.1 to 18.7
Units:	10 ⁹ /L	10 ⁹ /L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	2.9	2.9	0.0	S00015	3.5	3.6	0.1	S00029	7.4	7.8	0.4
S00002	5.3	5.3	0.0	S00016	2.7	2.9	0.2	S00030	6.8	6.9	0.1
S00003	5.2	5.6	0.4	S00017	5.0	4.8	-0.2	S00031	3.1	2.8	-0.3
S00004	4.7	5.2	0.5	S00018	5.3	5.3	0.0	S00032	14.7	14.7	0.0
S00005	2.7	2.6	-0.1	S00019	3.8	3.9	0.1	S00033	2.3	2.1	-0.2
S00006	5.0	5.2	0.2	S00020	5.8	6.3	0.5	S00034	5.4	6.6	1.2
S00007	2.4	2.5	0.1	S00021	4.6	4.9	0.3	S00035	2.6	2.6	0.0
S00008	3.2	3.5	0.3	S00022	4.0	4.3	0.3	S00036	16.7	18.7	2.0
S00009	5.4	5.4	0.0	S00023	3.4	4.3	0.9	S00037	2.4	2.5	0.1
S00010	2.9	3.0	0.1	S00024	4.5	4.0	-0.5	S00038	4.6	4.8	0.2
S00011	4.0	3.8	-0.2	S00025	3.6	3.8	0.2	S00039	2.3	2.4	0.1
S00012	6.0	5.8	-0.2	S00026	8.8	8.6	-0.2	S00040	4.1	4.3	0.2
S00013	5.1	5.4	0.3	S00027	4.5	4.4	-0.1				
S00014	4.9	4.7	-0.2	S00028	5.9	5.8	-0.1				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.