


## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	1.013 (0.996 to 1.030)	1.012 (0.995 to 1.028)
<b>Intercept:</b>	-0.007 (-0.027 to 0.013)	-0.006 (-0.025 to 0.014)
<b>Std Err Est:</b>	0.022	0.022

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R): 0.9988	Points (Plotted/Total): 40/40
Bias: 0.007	Outliers: None
XMean ± SD: 1.101 ± 0.434	Degrees Freedom: 38
YMean ± SD: 1.108 ± 0.439	Scatter Plot Bounds: None
Std Dev Diff: 0.023	
SubRange Bounds: None	

### Experiment Description

	X Method	Y Method
ExptDate:	11 Feb 2003	11 Feb 2003
Rep SD:	1	1
Result Ranges:	0.57 to 2.91	0.58 to 2.94
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_

Signature

Date

## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	0.99	1.02	0.03	S00015	1.00	0.97	-0.03	S00029	0.87	0.88	0.01
S00002	1.23	1.21	-0.02	S00016	0.99	1.02	0.03	S00030	0.72	0.73	0.01
S00003	0.60	0.60	0.00	S00017	0.81	0.83	0.02	S00031	1.18	1.18	0.00
S00004	2.31	2.34	0.03	S00018	1.19	1.15	-0.04	S00032	1.16	1.16	0.00
S00005	1.15	1.14	-0.01	S00019	1.06	1.09	0.03	S00033	1.04	1.08	0.04
S00006	1.09	1.08	-0.01	S00020	1.01	1.04	0.03	S00034	1.20	1.19	-0.01
S00007	0.98	0.96	-0.02	S00021	0.57	0.58	0.01	S00035	0.97	0.98	0.01
S00008	0.71	0.70	-0.01	S00022	0.95	0.98	0.03	S00036	1.25	1.26	0.01
S00009	1.08	1.06	-0.02	S00023	0.78	0.76	-0.02	S00037	1.29	1.32	0.03
S00010	1.05	1.07	0.02	S00024	1.03	1.04	0.01	S00038	0.68	0.68	0.00
S00011	0.89	0.91	0.02	S00025	2.91	2.94	0.03	S00039	0.82	0.87	0.05
S00012	2.06	2.10	0.04	S00026	1.07	1.07	0.00	S00040	1.09	1.12	0.03
S00013	1.15	1.12	-0.03	S00027	1.03	1.05	0.02				
S00014	1.23	1.21	-0.02	S00028	0.84	0.83	-0.01				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.