




Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	1.026 (0.994 to 1.059)	1.021 (0.989 to 1.054)
Intercept:	-3.6 (-12.2 to 4.9)	-2.4 (-11.0 to 6.2)
Std Err Est:	9.7	9.7

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9952	Points (Plotted/Total): 40/40
Bias: 2.8	Outliers: None
XMean ± SD: 243.6 ± 95.6	Degrees Freedom: 38
YMean ± SD: 246.4 ± 98.1	Scatter Plot Bounds: None
Std Dev Diff: 9.8	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	12 Feb 2003	12 Feb 2003
Rep SD:	1	1
Result Ranges:	13 to 644	13 to 656
Units:	10^9/L	10^9/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	170	162	-8	S00015	244	250	6	S00029	228	217	-11
S00002	360	359	-1	S00016	144	144	0	S00030	320	329	9
S00003	311	306	-5	S00017	194	195	1	S00031	229	230	1
S00004	223	228	5	S00018	125	132	7	S00032	644	656	12
S00005	214	211	-3	S00019	216	230	14	S00033	207	216	9
S00006	297	307	10	S00020	305	337	32	S00034	241	253	12
S00007	226	227	1	S00021	337	330	-7	S00035	173	173	0
S00008	316	310	-6	S00022	284	272	-12	S00036	13	13	0
S00009	289	286	-3	S00023	186	199	13	S00037	255	260	5
S00010	267	256	-11	S00024	270	269	-1	S00038	254	246	-8
S00011	142	139	-3	S00025	272	275	3	S00039	206	209	3
S00012	254	263	9	S00026	130	137	7	S00040	221	235	14
S00013	201	197	-4	S00027	213	212	-1				
S00014	375	405	30	S00028	189	182	-7				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.