


Alternate Method Comparison

X Method: BD

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	1.014 (0.997 to 1.030)	1.012 (0.996 to 1.029)
Intercept:	-0.038 (-0.113 to 0.037)	-0.032 (-0.107 to 0.043)
Std Err Est:	0.045	0.045

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9988	Points (Plotted/Total): 40/40
Bias: 0.023	Outliers: None
XMean ± SD: 4.434 ± 0.874	Degrees Freedom: 38
YMean ± SD: 4.457 ± 0.886	Scatter Plot Bounds: None
Std Dev Diff: 0.045	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	20 Dec 2002	20 Dec 2002
Rep SD:	1	1
Result Ranges:	2.55 to 7.91	2.55 to 8.01
Units:	x10 ¹² /L	x10 ¹² /L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	4.94	4.91	-0.03	S00015	4.15	4.17	0.02	S00029	4.29	4.31	0.02
S00002	3.94	3.99	0.05	S00016	5.15	5.16	0.01	S00030	3.51	3.50	-0.01
S00003	4.52	4.47	-0.05	S00017	3.76	3.78	0.02	S00031	6.24	6.27	0.03
S00004	4.47	4.44	-0.03	S00018	4.35	4.35	0.00	S00032	7.91	8.01	0.10
S00005	4.63	4.71	0.08	S00019	3.74	3.77	0.03	S00033	4.40	4.45	0.05
S00006	4.97	5.02	0.05	S00020	4.10	4.11	0.01	S00034	4.32	4.32	0.00
S00007	4.63	4.67	0.04	S00021	5.42	5.44	0.02	S00035	4.37	4.37	0.00
S00008	5.11	5.17	0.06	S00022	4.95	5.08	0.13	S00036	2.69	2.69	0.00
S00009	4.70	4.81	0.11	S00023	5.13	5.12	-0.01	S00037	4.22	4.23	0.01
S00010	4.19	4.23	0.04	S00024	2.55	2.55	0.00	S00038	4.22	4.31	0.09
S00011	4.28	4.28	0.00	S00025	4.15	4.16	0.01	S00039	4.56	4.6	0.04
S00012	4.62	4.56	-0.06	S00026	3.51	3.59	0.08	S00040	4.47	4.38	-0.09
S00013	4.44	4.44	0.00	S00027	4.07	4.12	0.05				
S00014	3.55	3.55	0.00	S00028	4.14	4.20	0.06				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.