


## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	1.014 (0.970 to 1.058)	1.005 (0.961 to 1.049)
<b>Intercept:</b>	0.0 (-0.5 to 0.4)	0.0 (-0.4 to 0.5)
<b>Std Err Est:</b>	0.5	0.5

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R): 0.9915	Points (Plotted/Total): 39/40
Bias: 0.1	Outliers: 1
XMean ± SD: 8.7 ± 4.0	Degrees Freedom: 37
YMean ± SD: 8.7 ± 4.1	Scatter Plot Bounds: None
Std Dev Diff: 0.5	
SubRange Bounds: None	

### Experiment Description

	X Method	Y Method
ExptDate:	08 Jan 2003	08 Jan 2003
Rep SD:	1	1
Result Ranges:	3 to 19	3 to 19
Units:	umol/L	umol/L
Analyst:	DIS	DIS
Comment:	SP#16 HEMOLYZED	

Accepted by: \_\_\_\_\_

Signature

\_\_\_\_\_

Date

## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias	
S00001	4	4	0	S00015	13	13	0	S00029	18	18	0	
S00002	5	5	0	S00016	O	9	2	-7	S00030	7	7	0
S00003	6	6	0	S00017	9	10	1	S00031	5	5	0	
S00004	9	9	0	S00018	8	8	0	S00032	9	9	0	
S00005	19	19	0	S00019	3	3	0	S00033	11	11	0	
S00006	7	7	0	S00020	6	6	0	S00034	9	9	0	
S00007	10	10	0	S00021	7	7	0	S00035	17	17	0	
S00008	12	12	0	S00022	6	7	1	S00036	4	4	0	
S00009	12	13	1	S00023	6	6	0	S00037	7	7	0	
S00010	7	7	0	S00024	7	5	-2	S00038	7	6	-1	
S00011	13	13	0	S00025	7	8	1	S00039	10	11	1	
S00012	5	5	0	S00026	6	7	1	S00040	17	17	0	
S00013	5	5	0	S00027	4	4	0					
S00014	11	11	0	S00028	10	10	0					

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.