




Alternate Method Comparison

X Method: BD

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	1.000 (0.927 to 1.074)	0.975 (0.902 to 1.049)
Intercept:	-0.1 (-5.5 to 5.2)	1.7 (-3.6 to 7.1)
Std Err Est:	1.0	1.0

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R):	0.9747	Points (Plotted/Total):	40/40
Bias:	-0.1	Outliers:	None
XMean ± SD:	72.7 ± 4.6	Degrees Freedom:	38
YMean ± SD:	72.5 ± 4.6	Scatter Plot Bounds:	None
Std Dev Diff:	1.0		
SubRange Bounds:	None		

Experiment Description

	X Method	Y Method
ExptDate:	11 Feb 2003	11 Feb 2003
Rep SD:	1	1
Result Ranges:	63 to 85	63 to 86
Units:	g/L	g/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	72	74	2	S00015	69	69	0	S00029	68	69	1
S00002	65	66	1	S00016	85	86	1	S00030	68	67	-1
S00003	68	67	-1	S00017	78	80	2	S00031	71	70	-1
S00004	75	74	-1	S00018	69	69	0	S00032	74	74	0
S00005	76	76	0	S00019	71	70	-1	S00033	75	75	0
S00006	73	74	1	S00020	70	70	0	S00034	75	75	0
S00007	75	75	0	S00021	69	68	-1	S00035	78	76	-2
S00008	72	72	0	S00022	70	70	0	S00036	81	81	0
S00009	74	73	-1	S00023	78	78	0	S00037	71	73	2
S00010	70	69	-1	S00024	75	74	-1	S00038	73	72	-1
S00011	69	69	0	S00025	69	69	0	S00039	69	70	1
S00012	68	69	1	S00026	76	73	-3	S00040	74	74	0
S00013	63	63	0	S00027	83	82	-1				
S00014	74	74	0	S00028	73	73	0				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.