


Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.995 (0.985 to 1.006)	0.995 (0.985 to 1.005)
Intercept:	-0.006 (-0.030 to 0.018)	-0.005 (-0.029 to 0.019)
Std Err Est:	0.034	0.034

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9995	Points (Plotted/Total): 40/40
Bias: -0.015	Outliers: None
XMean ± SD: 2.035 ± 1.071	Degrees Freedom: 38
YMean ± SD: 2.020 ± 1.066	Scatter Plot Bounds: None
Std Dev Diff: 0.034	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	11 Feb 2003	11 Feb 2003
Rep SD:	1	1
Result Ranges:	0.57 to 4.92	0.57 to 4.81
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	1.40	1.40	0.00	S00015	1.71	1.71	0.00	S00029	1.09	1.03	-0.06
S00002	4.92	4.81	-0.11	S00016	2.55	2.52	-0.03	S00030	1.73	1.67	-0.06
S00003	1.34	1.26	-0.08	S00017	2.42	2.49	0.07	S00031	1.13	1.11	-0.02
S00004	2.52	2.48	-0.04	S00018	3.64	3.65	0.01	S00032	4.29	4.31	0.02
S00005	1.71	1.68	-0.03	S00019	1.69	1.70	0.01	S00033	1.52	1.49	-0.03
S00006	1.85	1.83	-0.02	S00020	2.10	2.08	-0.02	S00034	0.82	0.82	0.00
S00007	1.30	1.29	-0.01	S00021	2.54	2.50	-0.04	S00035	2.59	2.51	-0.08
S00008	2.61	2.61	0.00	S00022	0.87	0.87	0.00	S00036	2.29	2.28	-0.01
S00009	1.27	1.26	-0.01	S00023	2.94	2.94	0.00	S00037	1.43	1.43	0.00
S00010	1.72	1.72	0.00	S00024	0.72	0.71	-0.01	S00038	4.60	4.59	-0.01
S00011	2.11	2.12	0.01	S00025	1.93	1.94	0.01	S00039	0.84	0.85	0.01
S00012	0.73	0.73	0.00	S00026	0.57	0.57	0.00	S00040	2.84	2.88	0.04
S00013	1.37	1.34	-0.03	S00027	3.25	3.17	-0.08				
S00014	1.15	1.15	0.00	S00028	3.30	3.29	-0.01				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.