




Alternate Method Comparison

X Method: BD

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	0.999 (0.986 to 1.012)	0.998 (0.985 to 1.011)
Intercept:	-1.8 (-6.9 to 3.2)	-1.6 (-6.6 to 3.5)
Std Err Est:	4.2	4.2

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9992	Points (Plotted/Total): 40/40
Bias: -2.3	Outliers: None
XMean ± SD: 381.4 ± 107.8	Degrees Freedom: 38
YMean ± SD: 379.1 ± 107.7	Scatter Plot Bounds: None
Std Dev Diff: 4.2	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	11 Feb 2003	11 Feb 2003
Rep SD:	1	1
Result Ranges:	172 to 681	175 to 680
Units:	umol/L	umol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	319	318	-1	S00015	510	512	2	S00029	323	321	-2
S00002	350	344	-6	S00016	393	391	-2	S00030	353	352	-1
S00003	361	356	-5	S00017	372	371	-1	S00031	376	373	-3
S00004	681	680	-1	S00018	318	314	-4	S00032	399	378	-21
S00005	322	315	-7	S00019	622	617	-5	S00033	217	217	0
S00006	429	429	0	S00020	280	280	0	S00034	335	332	-3
S00007	415	414	-1	S00021	382	375	-7	S00035	298	296	-2
S00008	436	435	-1	S00022	369	375	6	S00036	335	330	-5
S00009	213	208	-5	S00023	491	492	1	S00037	172	175	3
S00010	282	282	0	S00024	317	313	-4	S00038	356	352	-4
S00011	277	277	0	S00025	648	643	-5	S00039	354	354	0
S00012	419	423	4	S00026	382	379	-3	S00040	321	318	-3
S00013	416	414	-2	S00027	466	462	-4				
S00014	425	427	2	S00028	521	520	-1				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.