


## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	1.001 (0.967 to 1.035)	0.996 (0.961 to 1.030)
<b>Intercept:</b>	0.10 (-0.17 to 0.37)	0.14 (-0.13 to 0.41)
<b>Std Err Est:</b>	0.38	0.38

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R): 0.9946	Points (Plotted/Total): 40/40
Bias: 0.11	Outliers: None
XMean ± SD: 7.19 ± 3.57	Degrees Freedom: 38
YMean ± SD: 7.30 ± 3.57	Scatter Plot Bounds: None
Std Dev Diff: 0.37	
SubRange Bounds: None	

### Experiment Description

	X Method	Y Method
ExptDate:	20 Dec 2002	20 Dec 2002
Rep SD:	1	1
Result Ranges:	3.7 to 23.3	4.0 to 23.3
Units:	10 <sup>9</sup> /L	10 <sup>9</sup> /L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

## Alternate Method Comparison

**X Method: BD**

**Y Method: Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	4.8	4.8	0.0	S00015	5.0	5.1	0.1	S00029	9.6	10.1	0.5
S00002	8.9	8.8	-0.1	S00016	5.0	5.5	0.5	S00030	9.1	9.3	0.2
S00003	6.7	7.3	0.6	S00017	7.8	7.7	-0.1	S00031	6.8	6.3	-0.5
S00004	6.6	7.1	0.5	S00018	7.7	8.0	0.3	S00032	18.3	18.5	0.2
S00005	4.7	4.4	-0.3	S00019	5.1	5.1	0.0	S00033	4.2	4.0	-0.2
S00006	6.8	7.0	0.2	S00020	6.6	7.1	0.5	S00034	6.5	7.8	1.3
S00007	4.2	4.5	0.3	S00021	7.5	7.5	0.0	S00035	5.0	5.0	0.0
S00008	6.8	7.1	0.3	S00022	5.7	5.9	0.2	S00036	23.3	23.3	0.0
S00009	8.2	8.2	0.0	S00023	5.9	6.6	0.7	S00037	3.7	4.0	0.3
S00010	6.1	6.0	-0.1	S00024	6.0	5.2	-0.8	S00038	6.8	7.0	0.2
S00011	6.1	5.9	-0.2	S00025	5.1	5.3	0.2	S00039	4.1	4.3	0.2
S00012	9.5	9.1	-0.4	S00026	9.7	9.3	-0.4	S00040	6.9	7.0	0.1
S00013	8.0	8.4	0.4	S00027	6.1	6.0	-0.1				
S00014	5.8	5.6	-0.2	S00028	6.9	6.8	-0.1				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.