

making a difference



**FOR HYGIENIC
AND SAFE ROUTINE
URINE COLLECTION**

VACUETTE®

Urine System

The practical solution for
standardized preanalytics



greiner
BIO-ONE

A SYSTEMATIC SOLUTION

to cover routine needs.



OVERVIEW OF TUBE TYPES AND AREAS OF APPLICATION

VACUETTE® URINE TUBE WITHOUT ADDITIVE, ROUND BASE:

This tube can be used as a collection tube and/or transport container in the following areas of application:

- / Chemical urine tests*
- / Urine test strips*

It is recommended that the analysis is carried out as soon as possible after the collection. During the preanalytical phase, the stability of the sample is influenced by storage temperature and time. The cooler the storage temperature, the longer it can be stored. If the specimen cannot be tested within two hours of collection, it should be refrigerated (2-8°C) per CLSI GP16-A3 Guideline.

VACUETTE URINE TUBE WITHOUT ADDITIVE, CONICAL BASE:

The tubes are used as collection tube and/or transport container in the following area of application:

- / Urine microscopy*

The time between collection and analysis should not last longer than 2 hours. The recommended examination material for microscopical evaluation of sediment is mid-stream urine.



VACUETTE® URINE TUBES

round base (sterile)



VACUETTE® URINE CCM TUBES

conical base (sterile)



URINE CUP

with integrated transfer device (sterile or non-sterile)



URINE TRANSFER DEVICE

for the convenient sample transfer from the collection vessel



VACUETTE® URINE CCM TUBES WITH PRESERVATIVE (COUNT AND CULTURE MANNITOL):

The VACUETTE® Urine CCM Tubes contain a stabilizer in powder form, which keeps the bacterial count in the sample stable at room temperature (20-25°C) for up to 48 hours. The solubility of the additive stabilizes the urine sample quickly and is therefore espe-

cially suitable in cases when the analysis takes place later. This tube can be used as collection and/or transport container in the following area of application:

/ Microbiological tests

	Chemical urine analysis	Urine test strips	Microbiological studies	Urine sediment
Urine without preservative, round base	✓	✓	✓*	
Urine without preservative, conical base	✓	✓	✓*	✓
Urine CCM			✓	

* It is recommended that the analysis be carried out as soon as possible after collection.

THE OPTIMAL TUBE FOR MICROBIOLOGY

The VACUETTE® Urine CCM Tube keeps the bacterial count stable for up to 48 hours at room temperature so that there is no need to refrigerate.

The stabilizer in VACUETTE® Urine CCM Tubes (Count and Culture Mannitol Tube) is a combination of components including boric acid, sodium tetraborate, sodium formate and mannitol.

The quick solubility of the stabilizer is a distinguishing feature of the tube. In comparison with other additives, the improved solubility ensures immediate stabilization of the sample.

HOW SIGNIFICANT IS IMPROVED SAMPLE STABILITY IN LABORATORY ROUTINE?

- / Less re-collection
- / No need for refrigeration
- / Reliable diagnostics
- / Improved sample quality

shop.gbo.com

Find these and more products in our online shop

SAFETY CAP

completely leakproof and no splatter when opening

FILL-LEVEL MARK


precisely defined vacuum for reliable, correct urine to additive ratio

PET PLASTIC

virtually unbreakable specimen container

COUNT & CULTURE

mannitol powder



IMPROVED DIAGNOSTICS WITH THE VACUETTE® URINE CCM TUBE

Along with blood, urine is one of the most frequently analyzed specimen materials in the laboratory. Various analytical methods, such as microbiology culture, are used.

Some bacteria are very sensitive to changes in temperature or storage conditions. Therefore if the preanalytical conditions are not ideal, they can be difficult to detect.

A poor quality sample can result in incorrect analytical results and subsequent misdiagnosis. This requires the need to repeat analysis that, in addition to the impact to the patient, is associated with additional time re-

quirements and increased cost.

A common problem is bacterial growth in an unpreserved or unrefrigerated urine sample. If such an unpreserved urine sample is then stored or transported for extended periods of time, continued bacterial growth reduces the quality of the sample.¹

A urinary tract infection is typically defined as a bacterial count in urine $\geq 10^5$ CFU/ml (colony forming units per ml). A stable urine sample for microbiological analysis is important for an accurate diagnosis. If bacteriological urine examinations cannot be carried out in a timely manner (preferably within 2 hours^{2,3}), the sample has to be stored at 4°C or a preservative has to be added to act as a stabilizer.


The VACUETTE® Urine CCM Tube is ideal for samples that have long transportation and/or processing times.

RECOMMENDATIONS TO GUARANTEE A STABLE SAMPLE QUALITY:

- / In cases where the sample remains in the urine beaker for a longer period of time, it should be carefully and thoroughly mixed, thus homogenizing the sample prior to transferring.
- / Exclusive use of a sterile Urine Cup for urine collection will contribute to less risk of contamination, which could improve sample quality.
- / When using VACUETTE® Urine Tubes with stabilizers, tubes should be inverted several times so as to mix the urine and preservative homogeneously.
- / Recommended centrifugal conditions for conical base tubes: 5 min. / 400 g.

VACUETTE® Urine CCM Tubes reliably stabilize bacterial counts for up to 48 hours at room temperature.

- / Some bacteria divide every 20 minutes at room temperature.
- / If the sample quality is insufficient, pathogenic germs can be overgrown by contaminants.



PLEASE
OBSERVE THE
INSTRUCTIONS
FOR USE
at www.gbo.com

APPLICATION OF URINE CUP WITH INTEGRATED TRANSFER DEVICE*



1

If more than one tube will be filled, loosen the lid slightly more than a quarter turn to normalize the pressure inside the cup.



2

Peel back the safety label.



3

Insert the tube and wait until it is completely filled. Remove the tube afterwards.



4

Place the safety label in its original position and re-tighten the lid.

*For further collection devices, see table on page 11 and always refer to the correspondent Instructions for Use on www.gbo.com

ADVANTAGES



IMPROVED HYGIENE

Improved safety for laboratory and hospital personnel through the use of a closed system.



HIGHLY EFFICIENT AND COST SAVING

Collection, transportation, centrifugation, analysis and storage can all be executed in a single sterile tube.



SAFE AND CONVENIENT TRANSPORT

VACUETTE® Urine Tubes are hermetically sealed, sterile, virtually unbreakable and packed in compact racks.



HIGH QUALITY

Sterile VACUETTE® Urine Tubes have an exact dosage vacuum and meet highest quality standards.



CONTAMINATION RISK IS MINIMIZED

Contamination is nearly impossible when transferring into a secondary vessel.



PREANALYTICAL STANDARDIZATION

Errors produced as a result of collection, transportation and storage are reduced to a minimum.

PRODUCT OVERVIEW

VACUETTE® Urine Tubes

Item no.	Nominal volume	Cap colour	Ring colour	Tube size	Description	Packaging Inner	Packaging Outer
Z Urine No Additive - PREMIUM (Safety Twist Cap)							
456007	6 ml	● yellow	● yellow	13 x 100	Z Urine No Additive, Round Base	50 pcs.	1200 pcs.
Z Urine No Additive - Non-ridged (Pull Cap)							
455007	10 ml	● yellow	● yellow	16 x 100	Z Urine No Additive, Round Base	50 pcs.	1200 pcs.
455028	9 ml	● yellow	● yellow	16 x 100	Z Urine No Additive, Conical Base	50 pcs.	1200 pcs.
Urine CCM - Non-ridged (Pull Cap)							
454486	4 ml	● yellow	● black	13 x 75	Urine CCM, Round Base	50 pcs.	1200 pcs.
455243	9 ml	● yellow	● black	16 x 100	Urine CCM, Conical Base	50 pcs.	1200 pcs.

VACUETTE® Urine Accessories

Item no.	Name	Nominal volume	Description	Packaging Inner	Packaging Outer
Urine CCM Set (Urine CCM Tubes + Urine Transfer Device)					
453031	Urine CCM Set (13 x 75)	4 ml	VACUETTE® Urine CCM Tube + Urine Transfer Device (454486+450751), single-packed, Round Base	1 set	300 sets
453032	Urine CCM Set (16 x 100)	9 ml	VACUETTE® Urine CCM Tube + Urine Transfer Device (455243+450751), single-packed, Conical Base	1 set	300 sets
Urine Cup with Integrated Transfer Device					
724321	Urine Cup with Integrated Transfer Device 100 ml, yellow lid, screw cap, bulk-packed, sterile			-	300 pcs.
724322	Urine Cup with Integrated Transfer Device 100ml, yellow lid, screw cap, bulk-packed, non-sterile			-	300 pcs.
Urine Transfer Devices					
450251	Urine Transfer Device, short			50 pcs.	600 pcs.

REFERENCES:

1. VACUETTE® Preanalytics Manual p.49-53 Greiner Bio-One 2015
2. CLSI. Urinalysis; Approved Guideline - Third Edition. GP16-A3. Vol 29 No 4.
3. European Urinalysis Guidelines. Scan J Clin Lab Invest 2000. 60:1:1-96.

making a difference

www.gbo.com

GREINER BIO-ONE GMBH
KREMSMÜNSTER, AUSTRIA

PHONE +43 7583 6791-0
FAX +43 7583 6318
E-MAIL office@at.gbo.com



GREINER BIO-ONE IS A GLOBAL PLAYER.
FIND THE CONTACT DETAILS OF YOUR
LOCAL PARTNER ON OUR WEBSITE.



This product information is addressed exclusively to healthcare professionals. Devices of Greiner Bio-One are to be used by properly trained healthcare professionals only in accordance with the relevant Instructions for Use (IFU). For a listing of indications, contraindications, precautions and warnings, please refer to the Instructions for Use which accompanies each product or is available for download from our website at www.gbo.com (Download Center). For more information contact your local Greiner Bio-One sales representative or visit our website.

All information is provided without guarantee despite careful editing. Any liability, warranty or guarantee of Greiner Bio-One GmbH is excluded. All rights, errors and changes are reserved. If not stated otherwise, Greiner Bio-One GmbH has all copyrights and/or other (user-)rights in this/document / these documents, in particular to signs such as the mentioned (word-picture-)brands and logos. Any use, duplication or any other use of the rights of Greiner Bio-One GmbH is expressly prohibited.
Media owner: Greiner Bio-One GmbH / **Manufacturer:**

980005 en [rev. 14.04.2024]


greiner
BIO-ONE