



greiner bio-one

# Confirmation

## Transportation by aircraft

**Product description: Cryo.s™, Cryo.s™ with Datamatrix and Cryo.s™ Biobanking tubes**

**Item no.: 12x2xx, 12x2xx-xxx, 97x5xx, 97x5xx-xxx, 13x2xx, 13x2xx-xxx**

This Article meets the *pressure* requirements for *transportation by aircraft*

(Pressure requirements: packaging must be designed and constructed to prevent leakage that may be caused by changes in altitude and temperature during transportation aboard aircraft)

The hydrostatic pressure test was performed based on the ICAO.IATA DGR, 42nd Edition, Sections 6.3.5 and 5.0.2.9 and Title 49 CFR, Section 178.605


A sample of cryo vials was hydrostatically tested to 95 Kpa pressure differential as required for air shipment.

All tested cryo vials passed the test without any leakage

Certificate is valid until September 2016

21.09.2015

date

  
Quality Assurance

Greiner Bio-One GmbH  
Maybachstraße 2  
D-72636 Frickenhausen, Germany