



ENSURING SAMPLE INTEGRITY FROM COLLECTION TO ANALYSIS

TRANSPORT MONITORING

Precise Monitoring of Transport Bags or Samples

Tracie Transport Monitoring offers a modular solution for ensuring the integrity of biological samples during transport. Tracie ensures precise monitoring and data logging, crucial for maintaining sample quality and compliance with quality standards.

Key Features and Benefits

Transport Bag Tracking

Scan the transport bag once to capture essential tracking features, including location updates and temperature logging for the entire bag, ideal for straightforward transport needs

Sample Tracking

Each sample is scanned and assigned to the transport box, allowing for more detailed tracking, including real-time location tracking, individual temperature monitoring, and comprehensive data logging for each sample.

Temperature Control

Temperature data can either be monitored for the entire bag or for each sample, ensuring samples remain within required temperature ranges.

Integration with any existing System

Tracie Transport Monitoring can integrate with your LIS, allowing for smooth data transfer, but can also work as a standalone solution.

Compatible with any Transport Bag

Tracie works with any transport bag labeled with a unique barcode.




Tracie Transport Monitoring helps adhere to ISO 15189 and ISO 20658 standards by providing detailed documentation and traceability for every transport bag or each sample.



Tracie contributes to overall patient safety and reliable diagnostic outcomes by supporting a correct and secure transport of samples.



Tracie Transport Monitoring can be used as a standalone solution or in combination with Tracie Collect.



DETAILED MONITORING FOR ENHANCED TRANSPARENCY

Temperature Monitoring

Continuously track temperature of each bag or for each sample, with alerts for any deviations from specified temperature ranges.

Comprehensive Data Logging

Time stamps, Pickup/Checkout and Courier are logged at different stages during the process, providing a detailed record of the transport.


SIMPLIFYING LABORATORY PROCESSES

Efficient Check-In

Data transfer upon arrival of the transport bag in the lab.

Critical Data Display

View essential information such as time stamps and temperature curves, to help meet quality and compliance standards before processing.



GET STARTED WITH TRACIE TRANSPORT MONITORING!

Tracie Healthcare Solutions
www.tracie-health.com
service@tracie-health.com

Greiner Bio-One
www.gbo.com
office@gbo.com

LEARN MORE ABOUT TRANSPORT MONITORING

