




Erlenmeyer Shaker Flasks

Greiner Bio-One Erlenmeyer shaker flasks are ideal for suspension culture of mammalian cells, plants, microbes and microorganisms. The shaker flasks feature a patented 2-in-1 DuoCAP® which allows for flexibility between vented and non-vented applications. Each flask is manufactured utilizing an optically clear, inert polycarbonate resin which can be disposed of after single use or be made reusable when autoclaved.

<p>Product Specifications</p> <ul style="list-style-type: none"> ➤ 2-in-1 DuoCAP® for flexibility of sterile air exchange or tight seal ➤ Made of optically clear polycarbonate which is lighter and safer than glass ➤ Individually packaged and sterile to SAL 10⁻⁶ for long term culture applications ➤ Molded-in graduations for easy volumetric reference ➤ Non-pyrogenic, non-cytotoxic, and nuclease-free ➤ Available with flat or baffled base for increased agitation and aeration 	
---	--

Sterile Erlenmeyer Cell Culture Shaker Flasks					
Design	Cat. No.	Product Description		Pack Qty	Case Qty
Flat Bottom	679501	125mL		1	24
	679502	250mL		1	12
	679503	500mL		1	12
	679504	1000mL		1	6
	679505	2000mL		1	6
	679506	3000mL, Fernbach Style		1	3
Baffled Bottom	679511	125mL		1	24
	679512	250mL		1	12
	679513	500mL		1	12
	679514	1000mL		1	6
	679515	2000mL		1	6
	679516	3000mL, Fernbach Style		1	3