 DELTALAB	TECHNICAL DATA SHEET	309602E en	Date: 25/02/04
		Revision 3	Page: 1/2
GMDN code 10941			

Codes: 309602E, 409602
Screw cap container 25x90 mm with spoon

2003-2004 catalogue page: 146

Code	Case quantity	Case weight	Case volume	Cases per pallet
309602E	400	5.50	0,041	48
409602	400	4.90	0,041	48

Sanitary Product. IVDD 98 / 79 / CE

SPECIFIC USE AND DESCRIPTION:

Codes :	309602E	409602
Material :	CRYSTAL POLYSTYRENE	TRANSPARENT POLYPROPYLENE
Ext. diameter :	25.1 mm.	24.8 mm.
Ext. diameter with cap:	30.8 mm.	30.8 mm.
Total length (C) :	92.5 mm.	91.6 mm.
Maximum volume (D) :	31.7 ml. 29.7 ml with spoon.	31.1 ml. 29.1 ml with spoon.
Recommended volume (E) :	25 ml	25 ml
Graduations :	NO	NO
Withstand centrifugation :	7.000 g	10.000 g

- Conical bottom skirted container suitable for the transport of faeces and other solid samples from the collection site up to the laboratory where they will be analysed.
- Each container features a continuous thread and a blue ribbed polypropylene cap. Leakproof seal.
- Includes a removable spoon on underside of the cap so as to facilitate solid sample collection.
- The removable spoon is 6.9 cm long.
- Volume of the conical bottom: 1 ml.


Code 309602E :

- Container especially designed for hospital use.
- Labelled container. The label is designed to write the patient's name, the ward and its number, the type of specimen, the time and the date. Dimensions: 50 x 30 mm.

Elaborated by Marketing	Authorised by Management

 DELTALAB	TECHNICAL DATA SHEET	309602E en	Date: 25/02/04
GMDN code 10941		Revision 3	Page: 2/2

- The minimum ordering quantity: 30.000 units.



Elaborated by Marketing	Authorised by Management