

Certificate of Specification

Material : 84836.360 Alcool éthylique 96% GPR RECTAPUR

CHARACTERISTICS**SPECIFICATIONS**

Titre (V/V)	Min. 90,0 %
Dénatonium benzoate (11,8 -13,1 ppm)	Conforme
Méthyléthylcétone (V/V)	2,45 --> 2,99 %
Propanol-2 (V/V)	2,45 --> 2,99 %

We certify that this batch conforms to the specifications listed above.

BDL : Below detected limit.

Isabelle Habay, Head of laboratory - Briare
VWR International
Document printed on 02/2015

This document has been produced electronically and is valid without a signature.