

## NON STERILE STAINLESS STEEL HANDLES

Supplied in boxes of 10

### Handles in this range:

No.3 Graduated, 3L, No.7, No.9, B3, B3L

Fit all No.3 fitment blades

No.4 Graduated, 4L

Fit all No.4 fitment blades

| COMPONENT             | MATERIAL                                             |
|-----------------------|------------------------------------------------------|
| HANDLE                | Stainless Steel 17-4ph <b>OR</b> 430 Stainless Steel |
| UNIT PACK             | Cellophane Sleeve                                    |
| UNIT CONTAINER<br>BOX | Coated white lined chipboard                         |

All handles are manufactured in accordance with BS 2982 and BS EN 27740.

All Components are manufactured in the E.U.

All packets and boxes are marked with the current lot number and expiry dates for traceability.