

NON STERILE STAINLESS STEEL HANDLES

Supplied in boxes of 10

Handles in this range:

No.3 Graduated, 3L, No.7, No.9, B3, B3L

Fit all No.3 fitment blades

No.4 Graduated, 4L

Fit all No.4 fitment blades

COMPONENT	MATERIAL
HANDLE	Stainless Steel 17-4ph OR 430 Stainless Steel
UNIT PACK	Cellophane Sleeve
UNIT CONTAINER BOX	Coated white lined chipboard

All handles are manufactured in accordance with BS 2982 and BS EN 27740.

All Components are manufactured in the E.U.

All packets and boxes are marked with the current lot number and expiry dates for traceability.