

3 PART SYRINGES- Luer Slip

Reference Items	Description
REF:	SY3-XSC-EC
Associated group	Injection Systems
Product description	Syringe
Intended use	Injection therapy



Reference code	Description	Size (ml)	Type	Scaling (ml)
SY3-2SC-EC	3 part syringe 2 ml luer slip	2	CENTRIC	0.1
SY3-5SC-EC	3 part syringe 5 ml luer slip	5	CENTRIC	0.2
SY3-10SC-EC	3 part syringe 10 ml luer slip	10	CENTRIC	0.2
SY3- 10ESC-EC	3 part syringe 10 ml luer slip	10	ECCENTRIC	0.2
SY3-20ESC-EC	3 part syringe 20 ml luer slip	20	ECCENTRIC	1
SY3-30SC-EC	3 part syringe 30 ml luer slip	30	ECCENTRIC	1



2- 5- 10- 20 ml



30 ml





PRODUCT DATA SHEET

<u>Product Registrations</u>	
CE- marked	CE
Classification	Medical Devices
Medical Product class	Class I
Rule	Rule 2
Guideline	MDD 93/ 42/ EEC
UMDNS	13929
<u>Standard norms:</u>	
ISO 7886-1/ EN ISO 15223-1/ ISO 11607-1/ ISO 11607-2/ EN ISO 9001 EN ISO 13485/ ISO 10993-1/ ISO 10993-4/ ISO 10993-5/ ISO 10993-7/ EN 1707/ ISO 594-1/ ISO 594-2	


<u>Manufacturing details</u>	
Legal Manufacturer	Nipro Europe
Country of origin	China

<u>Product materials :</u>	
All materials used to produce our 3 part syringes are medical grade, latex-Free and PVC-Free.	
Poly-propylene	Plunger, Barrel
Isoprene rubber	Gasket
Silicone	Lubricant


<u>Product Sterilization</u>		
Sterile	ETO Sterilization procedure Ethylene – Oxide	
Do NOT re-use	Single use only	
<i>Each outer box (or shipping carton) contains a chemical indicator, which ensures sterility. This indicator is a blue sticker, which turns red when sterilized.</i>		

PRODUCT DATA SHEET

Packaging

- LOT-number
- Expiry : 5 years 
- REF: SY3-XSC-EC
- QTY- pieces
- Barcode

Example of Label:



An individual package, is sterile paper for blister package. Inner box is made of duplex board, and outer box of double corrugated board.

Transport conditions

Closed and dry

Storage conditions

Only open the packaging before use to guarantee sterility.

Reference	Size (ml)	Scaling	Quantity inner box	Quantity Outer box	Net Weight	Gross Weight	INNER BOX			OUTER BOX			
							L mm	W mm	H mm	L mm	W mm	H mm	M ³
SY3-2SC-EC	2	0.1	100	3200	5,000	16,000	420	125	70	580	440	530	0,135
SY3-5SC-EC	5	0.2	100	2400	6,660	16,000	420	125	93	580	440	530	0,135
SY3-10SC-EC	10	0.2	100	1500	9,330	14,000	420	165	110	580	440	530	0,135
SY3-10ESC-EC	10	0.2	100	1500	9,830	45,000	420	165	110	580	440	530	0,135
SY3-20SC-EC	20	1	50	800	15,000	12,000	215	210	140	580	440	460	0,101
SY3-30SC-EC	30	1	50	600	9,830	5,900	285	180	140	580	380	460	0,101

15 Languages on packaging according the specifications of European Medical Device Directive:

- English, French, Dutch, German, Spanish, Italian, Portuguese, Greek, Swedish, Danish, Norwegian, Finnish, Russian, Turkish, Polish



Contact details

Nipro Europe NV
Weihoek 3H
1930 Zaventem
Belgium

Phone: +32 (0)2 725 55 33 / Fax: +32 (0)2 725 70 41

E-mail: info@nipro-mediceurope.com

Web: www.nipro-mediceurope.com