



## Product Specification Sheet

---

### VWR RACK 80-Well, PP

---

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

*VWR International certifie que les produits mentionés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

| Cat. No. | For tubes  | No. of holes | Colour   | WxDxH        | Pk |
|----------|------------|--------------|----------|--------------|----|
| 211-0204 | 1.5&2.0 ml | 80           | Natural  | 225x67x28 mm | 5  |
| 211-0205 | 1.5&2.0 ml | 80           | Orange   | 225x67x28 mm | 5  |
| 211-0206 | 1.5&2.0 ml | 80           | Yellow   | 225x67x28 mm | 5  |
| 211-0207 | 1.5&2.0 ml | 80           | Blue     | 225x67x28 mm | 5  |
| 211-0208 | 1.5&2.0 ml | 80           | Green    | 225x67x28 mm | 5  |
| 211-0209 | 1.5&2.0 ml | 80           | Pink     | 225x67x28 mm | 5  |
| 211-0210 | 1.5&2.0 ml | 80           | Assorted | 225x67x28 mm | 5  |

**Material:** Polypropylene

**Properties:** Good chemical resistance

**Temperature Resistance:** 0° up to +121°C

**Feature:** Autoclaveable

Injection-molded of polypropylene homopolymer, which are inert or minimally reactive with most common laboratory reagents. PP used conforms to current REACH requirements.

**Disclaimer:** VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us. Date of issue: **May 2014**