



Certificate of Conformity

ISOLATION GOWN 60G - 115X136 CM – WHITE 113-8321

CHARACTERISTICS:

The disposable gown with back fastening is made of non woven breathable SMS fabric. Ideal for use in healthcare environments and hospitals.

- Long sleeves, knitted cuffs
- Round neck
- Ties on neck and waist

PHYSICAL PROPERTIES

Material	Non-woven SMS 60grams/m2	
Size	One size 115*136cm	
Length:	115cm +/- 2cm	
Width:	136cm +/- 2cm	
Sleeve length:	57cm +/- 2cm	
Pockets:	no pockets	
Colour	White	
Tensile Strength	MD: >60N/5cm	CD: >45N/5cm
Ultimate Elongation	>30 %	

PACKING:

1pcs/bag,50bags/carton

QUALITY ASSURANCE CRITERIA:

SMS gowns are guaranteed to meet 4.0 AQL level II sampling plan per MIL-std-150E for the following defect criteria:

- | | |
|------------------------------|--|
| 1. Foul odor | 4. Holes/tears in fabric material caused by sewing machine |
| 2. Irregular fabric material | 5. Physical dimensions non-compliance |
| 3. Dirty fabric material | |



Certificate of Conformity

***Disclaimer:** VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.*

Date of issue: **Feb 2020**