

| | | | |
|---|-----------------------------|--------------------------|---------------------|
|  | TECHNICAL DATA SHEET | M 5555B | Date: 24/03/14 |
| | | Revision 0 | Page: 1/3 |

**Codes: M 5555B, M 6080B, M 80120B, M 70110B, M 100150B,
M 120180, M 160220B, M 180250B, M 200300B**

White strip zip-lock bags

2014 Catalogue Page : 112

| Code | Gauge | Dimensions Mm | Case Quantity | Case Weight | Case Volume |
|------------------|-------|------------------|------------------|----------------|----------------|
| M 5555B | 200 | 55 x 55 | 2,000 | 1.13 | 0.0028 |
| M 6080B | 200 | 60 x 80 | 1,000 | 0.67 | 0.0026 |
| M 80120B | 200 | 80 x 120 | 1,000 | 1.27 | 0.0037 |
| M 70110B | 200 | 70 x 110 | 1,000 | 0.95 | 0.0037 |
| M 100150B | 200 | 100 x 150 | 1,000 | 1.59 | 0.0044 |
| M 120180B | 200 | 120 x 180 | 1,000 | 2.35 | 0.0065 |
| M 160220B | 200 | 160 x 220 | 1,000 | 4.32 | 0.0080 |
| M 180250B | 200 | 180 x 250 | 1,000 | 4.43 | 0.0125 |
| M 200300B | 200 | 200 x 300 | 1,000 | 6.82 | 0.0145 |

DESCRIPTION

Made of transparent polyethylene, very resistant to shocks and tearing. For all solid samples.

Suitable for laboratories (for tubes, swabs, etc.) as well as industries (for spare pieces, small electronic components, buttons, etc.).

Features white strip for writing and better sample identification.



| | | |
|---|---|---|
| Written | Verified | Authorised |
|  |  |  |