

SPECIAL
DISCOUNT*


greiner
BIO-ONE

**CLEAR INSIGHTS
FOR BETTER
OUTCOMES**

μCLEAR® & UV-STAR® PLATES

Clear-bottom microplates for optical readouts

Experience Precision with **μClear®** Plates!

Discover exceptional clarity and reliability in optical readouts with **μClear®** plates. Featuring a 190 μm film bottom, **μClear®** plates offer high spectral transmission, making them perfect for reader-based applications and standard microscopic applications at low to medium magnification. Benefit from low autofluorescence, a high signal-to-noise ratio, and high bottom planarity, ensuring reliable and reproducible results every time.

For measurements in the lower UV range, such as the determination of nucleic acid and protein concentrations at 260 nm or 280 nm, **UV-Star® film bottom microplates** are essential. Made entirely from cycloolefin, these plates provide excellent transmission down to 230 nm.

KEY FACTS:

- / 96 and 384 well format for bottom reading
- / High image quality
- / Low autofluorescence
- / High signal-to-noise ratio
- / Low autofocus times
- / Transmission in the UV with UV-Star® plates
- / Surfaces for biochemical and cell-based assays



REQUEST YOUR INDIVIDUAL OFFER
OR SAMPLE FOR TESTING



Greiner Bio-One GmbH Frickenhausen, Germany
PHONE +49 7022 948-0 / FAX +49 7022 948-514 / E-MAIL info@de.gbo.com
Greiner Bio-One is a global player.
Find the contact details of your local partner on our website.

www.gbo.com

PRODUCT OVERVIEW

REQUEST YOUR INDIVIDUAL OFFER
OR SAMPLE FOR TESTING

μCLEAR® FILM BOTTOM PLATES

Bottom: μClear® film bottom, Material: PS (polystyrene)

Item no.	Well format	Well profile	Product colour	Surface treatment	Lid	Sterile	Qty. inner / outer
655090	96	F-bottom / Chimney Well	● black	TC	yes, condensation rings	+	8 / 32
655095	96	F-bottom / Chimney Well	○ white	untreated	no		10 / 40
655096	96	F-bottom / Chimney Well	● black	untreated	no		10 / 40
655098	96	F-bottom / Chimney Well	○ white	TC	yes, condensation rings	+	8 / 32
781090	384	F-bottom	● black	TC	yes	+	20 / 120
781091	384	F-bottom	● black	TC	yes	+	8 / 32
781098	384	F-bottom	○ white	TC	yes	+	8 / 32

Further μClear® plates are available on request.

UV-Star® MICROPLATES

Bottom: COC (Cycloolefin Copolymer) film, Product colour: clear,
Material: COC, Surface treatment: untreated, Lid: no

Item no.	Well format	Well profile	Sterile	Qty. inner / outer
655801	96	F-bottom / Chimney Well		10 / 40
781801	384	F-bottom		10 / 40

Further UV-Star® plates are available on request.



GET YOUR INDIVIDUAL DISCOUNT AND/OR SAMPLE TODAY!*

* This discount is available for new customers on their first order of this product in specific regions.

For inquiries, please use the contact form accessible via the QR code or link.